

Programming DBS: Trouble Shooting Strategies

Susan Heath, M.S., R.N.
Movement Disorders Clinical Nurse Specialist
Associate Director of Clinical Care
San Francisco VA
Parkinson's Disease Research and Clinical Center
(PADRECC)



Common Problems Encountered in DBS Therapy

Lead location → suboptimal placement?

Programming suboptimal?

Device-related issues or malfunction?

Has there been an over zealous medication
reduction?

(pt under-medicated?) Or unrealistic expectations?


Is it a Lead Location Problem?



- What does the post-op MRI show?
- 1.5 tesla magnet only

Paul S. Larson, R. Mark Richardson, Philip A. Starr, Alastair J. Martin (2008).

Magnetic Resonance Imaging of Implanted Deep Brain Stimulators: Experience in a Large Series. Stereotactic Functional Neurosurgery; 86:92–100

Post-Operative MRI





Lead Location Directly Influences Programming: STN






FIGURE 2



Starr, Theodosopoulos, and Turner.
Surgery of the subthalamic nucleus: use of movement-related neuronal activity for surgical navigation.
Neurosurgery. 2003 Nov;53(5):1146-9; discussion 1149



Common Problems Encountered in DBS Therapy

Lead location → suboptimal placement?

Programming suboptimal?

Device-related issues or malfunction?

Has there been an over zealous medication reduction?
(pt under-medicated?) Or unrealistic expectations?

Programming with Gusto

- Has each contact been screened using *maximum tolerated voltage* to assess ceiling threshold?
 - Determines therapeutic window
 - Titrate just below maximum tolerated voltage
- Kinetra batteries require higher voltages



National VA Parkinson's Disease
CONSORTIUM
Education • Collaboration • Advocacy

Three Dimensional STN Anatomy

- Medial and posterior = lemniscal fibers
- Lateral and anterior = capsular fibers
- Anterior to motor = limbic pathways ?
- Deep and medial = SNr = limbic
- Medial and deep = IIIrd nerve
- Lateral to STN, genu of capsule = conjugate gaze eye signs



National VA Parkinson's Disease
CONSORTIUM
Education • Collaboration • Advocacy

STN Side Effects at Low Thresholds

Expected side effects from lead location:

- Persistent Paresthesias → Too posterior or medial (Lemniscal fibers)
- Dysarthria, pulling/contractions → Too lateral (Capsule)
- Diplopia → Too deep, **anteromedial**
- Dyskinesias → In optimal location
- Mood changes → Too anterior and medial
- **No effect at high voltages** → Too superior or anterior

Common Problems Encountered in DBS Therapy

Lead location → suboptimal placement?

Programming suboptimal?

Device-related issues or malfunction?

Has there been an over zealous medication reduction?

(pt under-medicated?)

Is Problem a Device Malfunction?

- Impedance checks – verify electrical integrity
- Kinetra (use 4.0V or 4.5V to run checks)
- Soletra (use std 1.5V to run checks)
- Short circuit = two contacts or wires connecting → rapid battery depletion
- Open Circuit = loss of one or more contacts with sudden loss of efficacy

Impedance Readings

– Impedance & current drain indicate status of system (applies to Soletra)

- Normal: impedance 600-1300 Ω ; current drain <12-30 μA
- Open Circuit:
 - Soletra impedance >2000 Ω ; current drain <7-9 μA
 - Kinetra impedance >4000 Ω ; current drain <7-9 μA
- Short Circuit: impedance <250 Ω ; current drain >500 μA

Open Circuit - loss of one or more contacts

Electrode Impedances
(1.5 V, 210 μ s, 30 Hz)

Electrode Pair	Impedance (Ω)	Current (μ A)
0 & c	>2000	<7
1 & c	>2000	<7
2 & c	>2000	<7
3 & c	>2000	<7
0 & 1	1137	14
0 & 2	1480	12
0 & 3	1535	12
1 & 2	1159	14
1 & 3	1381	13
2 & 3	1217	13

Wire Bond failure =
Loss of unipolar, but
preserved bipolar
function.

This pt
had loss of
efficacy in
Monopolar
setting, but,
preserved function
in bipolar setting

Lead Fracture – Open Circuit



National VA Parkinson's Disease
CONSORTIUM
Education • Collaboration • Advocacy

Short Circuit

- Short Circuit: Two wires making contact
- Impedance $<250\ \Omega$; current drain $>500\ \mu\text{A}$
- Rapid depletion of battery life



National VA Parkinson's Disease
CONSORTIUM
Education • Collaboration • Advocacy


Kinetra Impedance Challenges

- Less reliable and harder to interpret
- Run 900's -1800's (wider range, and sometimes higher than Soletra)
- Often similar readings between many contacts/settings – unclear how to interpret



National VA Parkinson's Disease
CONSORTIUM
Education • Collaboration • Advocacy

Kinetra Battery with 'fluid in the boot'



National VA Parkinson's Disease
CONSORTIUM
Education • Collaboration • Advocacy

Low Battery Issues

- **Soletra** low battery = $< 3.60V$, common for pts to want 'more juice' or feel less well
- **Kinetra** low battery = $< 2.40V$, uncommon for pts to feel loss of efficacy (start surgical referral at 2.50V)
- Advanced PD patients have required inpatient admission due to severe PD sx's after battery failed
- Goal is to monitor battery and replace *before* battery fails
- Need 2-4 weeks to plan replacement surgery



National VA Parkinson's Disease
CONSORTIUM
Education • Collaboration • Advocacy

Common Problems Encountered in DBS Therapy

Lead location → suboptimal placement?

Programming suboptimal?

Volts pushed up to just below side effects?

At high voltages, what are the side effects?

Device-related issues or malfunction?

- Hardware issues: wire fracture or crushed?
- Factory recall :‘wire bond’ disruption in battery

Has there been an over zealous medication reduction?
(pt under-medicated?)

Patient Education: Identify Secret Unrealistic Expectation(s)

- Managing Pt expectations “mantra”
 - DBS is not a cure for PD, can’t find ‘holy grail’ settings
 - Optimal results from stimulation different for each patient
- DBS settings and medications often need to be adjusted concurrently with STN (goes against rule to change one variable at a time)

Adjunctive therapies needed for advanced pts: physical therapy, occupational therapy and speech therapy often indicated

- Depression common after medication reduction



National VA Parkinson's Disease
CONSORTIUM
Education • Collaboration • Advocacy